


# β-Shielding

for use with Modular-Lab PharmTracer

## β-shielding devices

A wide range of β-shielding accessories (Plexiglas®) is available for Modular-Lab PharmTracer increasing radiation safety when working with β-emitters.



### Transport box and upper shell

The waste box (transport box) ensures the safe storage of the contaminated cassettes during the decay of the activity. The upper shell is attached to the waste box to guarantee a safe cassette-drop into the transport box.

Waste Box	
Shielding material:	Plexiglas®
Dimensions (L x W x H):	50.7 x 44.7 x 21.9 cm
Weight:	12.5 kg
Article No.:	3111-0666

Upper Shell	
Shielding material:	Plexiglas®
Dimensions (L x W x H):	50.7 x 44.7 x 17.6 cm
Weight:	4.2 kg
Article No.:	3111-0665

### β-shielding for <sup>90</sup>Y and <sup>177</sup>Lu cassettes

The cassette shielding for disposable <sup>90</sup>Y and <sup>177</sup>Lu cassettes is attached to your Modular-Lab PharmTracer. It minimizes the user's contact with the contaminated cassette during disposal.

Shielding material:	Plexiglas®
Dimensions (L x W x H):	46.9 x 17.7 x 39.7 cm
Weight:	2.7 kg
Article No.:	3111-0663



### Vial downholder

The downholder facilitates the insertion of a needle into the vial.

Dimensions (L x W x H):	11.3 x 9.3 x 10.4 cm
Weight:	0.6 kg
Article No.:	3111-0673

# β-Shielding



## β-shielding container for product vial

For safe transport of a 30 ml product vial

Shielding material:	Plexiglas®
Dimensions (L x W x H):	12.9 x 8.8 x 14.5 cm
Weight:	1.1 kg
Article No.:	3111-0671

## β-shielding container for waste vial 100 ml

For safe transport of a 100 ml waste vial

Shielding material:	Plexiglas®
Dimensions (L x W x H):	15.1 x 11 x 14.5 cm
Weight:	1.3 kg
Article No.:	3111-0672

## β-shielding gripper for product vial

For safe handling of a product vial

Shielding material:	Plexiglas®
Dimensions (L x W x H):	28.5 x 13 cm
Weight:	0.3 kg
Article No.:	3111-0674

## β-shielding gripper for waste vial

For safe handling of a waste vial

Shielding material:	Plexiglas®
Dimensions (L x W x H):	28.5 x 13 cm
Weight:	0.3 kg
Article No.:	3111-0675

The accessories presented herein are suitable for β-shielding only!

### Eckert & Ziegler Eurotope GmbH

Robert-Rössle-Strasse 10  
13125 Berlin  
Germany  
Phone: +49 30 941084 197  
Fax: +49 30 941084 470

eurotope@ezag.de  
www.radiopharma.com

### Your contact in USA & Canada: Eckert & Ziegler Radiopharma, Inc.

25 Upton Drive  
Wilmington, MA 01887  
USA  
Phone: +1 508 497 0060  
Fax: +1 508 497 0061

infoRAU@ezag.com  
www.radiopharma.com