

Modular-Lab eazy

Cassette-based solution for routine production



Cassette Diversity

The fully automated synthesis system Modular-Lab eazy works with sterile, disposable cassettes.

The use of sterile, disposable synthesis cassettes ensures the easy handling through click'n'start technology. Just by changing the cassette the Modular-Lab eazy is ready for the next synthesis run - even if it is a different tracer or nuclide than the one the customer has labeled just before. Due to the cassettes' one-time use no reconditioning and cleaning or sanitation routine are necessary.

Cassettes are assembled under GMP-compliant clean room conditions and gamma-sterilized. All cassette parts were tested for their suitability for the specific syntheses. All materials are free of animal derived ingredients. A shelf life of 18 months can be guaranteed.

Cassettes are available for the following ^{68}Ga application:

	Cassette for Modular-Lab eazy	Article No.
	Cassette for ^{68}Ga Multi-Generator-Elution	C1-GA68-MULTI
	Tubing set for elution of registered $^{68}\text{Ge}/^{68}\text{Ga}$ generator for GalliaPharm®	ML-EAZY-GALLIA-TUBES
^{68}Ga	^{68}Ga DOTA-peptides (acetone free pre-purification) with CM cartridge post-purification	C0-GA-PEP
	^{68}Ga labeled PSMA ^{HBED} with CM cartridge post-purification	C0-GA-PSMA11
	^{68}Ga labeled PentixaFor with SCX cartridge pre-purification and CM cartridge post-purification	C0-GA-PXF
	^{68}Ga labeled FAPI-46 with SCX cartridge pre-purification and CM cartridge post-purification	C0-GA-FAP
^{90}Y	^{90}Y labeled PSMA-617 with CM cartridge post-purification	C0-YDOTAPSMA-CM
	^{90}Y DOTA-peptides with CM cartridge post-purification	C0-YDOTAPEP-CM
^{177}Lu	^{177}Lu labeled PSMA-617 with CM cartridge post-purification	C0-LUDOTAPSMA-CM
	^{177}Lu DOTA-peptides with CM cartridge post-purification	C0-LUDOTAPEP-CM

**Eckert & Ziegler
Eurotope GmbH**

Robert-Rössle-Strasse 10
13125 Berlin
Germany
Phone: +49 30 941084 197
Fax: +49 30 941084 470

eurotope@ezag.de
www.radiopharma.com

Your contact in USA & Canada:
Eckert & Ziegler Radiopharma, Inc.

25 Upton Drive
Wilmington, MA 01887
USA
Phone: +1 508 497 0060
Fax: +1 508 497 0061

infoRAU@ezag.com
www.radiopharma.com

 **Eckert & Ziegler**