

Rack for stacking two $^{68}\text{Ge}/^{68}\text{Ga}$ Generators



DoubleGen

Save valuable space while working with two $^{68}\text{Ge}/^{68}\text{Ga}$ generators

Safe and compact arrangement of two generators

DoubleGen allows you to safely stack two of our Eckert & Ziegler $^{68}\text{Ge}/^{68}\text{Ga}$ generators on top of each other. This way the footprint used up by your generators is dramatically reduced to open up room in your work environment.

The two generators are mounted separately on the DoubleGen, so they can be removed independently of one another.

Safe storage of the upper generator is guaranteed through fixing of the handles with a dedicated clamp allowing for easy installation and removal of the generator.



Convenient storage of your HCl bags

For further efficiency, the DoubleGen enables you to install the required hydrochloric acid bags directly above your generators.

Technical Data

Dimensions	150 x 200 x 680 mm (W x D x H) 5.9" x 7.9" x 26.8" (W x D x H)
------------	---

Including one GalliaPharm® generator

Weight	19 kg / 41.9 lbs
--------	------------------

Point load	185 N
------------	-------

Surface load	0.73 N / cm ²
--------------	--------------------------

Including two GalliaPharm® generators

Weight	33 kg / 72.8 lbs
--------	------------------

Point load	325 N
------------	-------

Surface load	1.3 N / cm ²
--------------	-------------------------



Scan the QR code and discover all our generator accessories!



Eckert & Ziegler Eurotope GmbH

Robert-Rössle-Strasse 10
13125 Berlin
Germany
Phone: +49 30 941084 197
Fax: +49 30 941084 470

eurotope@ezag.de
www.radiopharma.com

Your contact in USA & Canada: Eckert & Ziegler Radiopharma, Inc.

25 Upton Drive
Wilmington MA 01887
USA
Phone: + 1 508 497 0060
Fax: + 1 508 497 0061

InfoRAU@ezag.com
www.radiopharma.com